



# Welch Allyn FlexiPort™ Disposable Blood Pressure Cuffs

True Standardisation in a Snap.

**WelchAllyn®**

Advancing Frontline Care™

# Simple to order. Simple to use. Simply amazing.

## Welch Allyn FlexiPort™ Disposable Blood Pressure Cuffs

True standardisation in a snap.

Convenient, comfortable, and affordable, Welch Allyn disposable blood pressure cuffs are the top choice for infection-control-conscious hospitals and clinics around the world. That's because we offer virtually every type of cuff an institution may need, along with the ability to simplify inventory like no other cuff provider with our patented FlexiPort connection.

With FlexiPort Disposable Cuffs, it's easy to move patients from room to room, floor to floor, and department to department. Simply equip all devices with FlexiPort fittings, assign one cuff to your patients as they enter the hospital, and connect the patient's cuff to the device.

### When you standardise on Welch Allyn FlexiPort Disposable Blood Pressure Cuffs:

- Every Cuff can be used as one- or two-tube cuff
- Part numbers are reduced by up to 60 percent, simplifying ordering and inventory
- You get every cuff size needed—from the tiniest neonatal\* through thigh cuffs
- You reduce the risk of cross-contamination

\*Neonate cuffs are not FlexiPort compatible



Patented FlexiPort™ technology allows one cuff to work in virtually any patient care area.

### The Disposable Difference

By converting to disposable cuffs, you can significantly reduce the risk of cross-contamination and related infections.

Reusable blood pressure cuffs can carry such bacteria as *Klebsiella pneumoniae*, Norwegian Scabies, *E. cloacae*, and *Staphylococcus aureus*—bacteria potentially deadly to infants and adults alike.

Reused blood pressure cuffs are a primary carrier of MRSA, which is responsible for 12% of all hospital-related infections in the U.S.



## What's FlexiPort?

The FlexiPort fitting provides a single point of connection between devices and cuffs, making any cuff a one- or two-tube cuff.



Rotatable port reduces stress to cuff tubing and port, improving patient comfort and cuff durability



Colour coded for fast, accurate selection of the correct cuff size

Meets all the latest clinical guidelines for proper fit from AAMI and AHA

No sewn seams to tear or fray  
Soft, latex-free material minimises risk of allergic reactions and ensures patient comfort

Unique ability to connect and disconnect tubing making one- and two-tube cuffs interchangeable

# Reduce Part Numbers by 60%



Port-2

Port-1

## The Only True Standardisation Solution

With FlexiPort, Welch Allyn is the only company to offer a true cuff standardisation solution for your facility. By connecting tubes at the cuff, disposable FlexiPort cuffs not only make it easier for caregivers to find and change the right cuff, but also reduce part numbers by up to 60 percent, simplifying inventory and ordering.

Contact your local Welch Allyn sales representative for more information on simplifying your cuff inventory with Welch Allyn.

		Option 1	Option 2							
		Attach above fittings directly to your current device tubing.	Buy FlexiPort cuffs with tubing and connectors to fit your blood pressure devices without making changes.							
		No Tubes or Connectors Required	One-Tube				Two-Tube			
Size (20 cuffs/case)	Range	FlexiPort	Screw	Tri-Purpose	Locking	Bayonet	Screw	Locking	Tri-Purpose & Empty Tube	Inflation Sys.
Small Infant (6)	7-10 cm	SOFT-06	SOFT-06-1SC	SOFT-06-1TPE	SOFT-06-1MQ	SOFT-06-1HP	n/a	n/a	n/a	n/a
Infant (7)	9-13 cm	SOFT-07	SOFT-07-1SC	SOFT-07-1TPE	SOFT-07-1MQ	SOFT-07-1HP	SOFT-07-2SC	SOFT-07-2MQ	SOFT-07-2TPE	SOFT-07-2BV
Small Child (8)	12-16 cm	SOFT-08	SOFT-08-1SC	SOFT-08-1TPE	SOFT-08-1MQ	SOFT-08-1HP	SOFT-08-2SC	SOFT-08-2MQ	SOFT-08-2TPE	SOFT-08-2BV
Child (9)	15-21 cm	SOFT-09	SOFT-09-1SC	SOFT-09-1TPE	SOFT-09-1MQ	SOFT-09-1HP	SOFT-09-2SC	SOFT-09-2MQ	SOFT-09-2TPE	SOFT-09-2BV
Small Adult (10)	20-26 cm	SOFT-10	SOFT-10-1SC	SOFT-10-1TPE	SOFT-10-1MQ	SOFT-10-1HP	SOFT-10-2SC	SOFT-10-2MQ	SOFT-10-2TPE	SOFT-10-2BV
Adult (11)	25-34 cm	SOFT-11	SOFT-11-1SC	SOFT-11-1TPE	SOFT-11-1MQ	SOFT-11-1HP	SOFT-11-2SC	SOFT-11-2MQ	SOFT-11-2TPE	SOFT-11-2BV
Adult Long (11L)	25-34 cm	SOFT-11L	SOFT-11L-1SC	SOFT-11L-1TPE	SOFT-11L-1MQ	SOFT-11L-1HP	SOFT-11L-2SC	SOFT-11L-2MQ	SOFT-11L-2TPE	SOFT-11L-2BV
Large Adult (12)	32-43 cm	SOFT-12	SOFT-12-1SC	SOFT-12-1TPE	SOFT-12-1MQ	SOFT-12-1HP	SOFT-12-2SC	SOFT-12-2MQ	SOFT-12-2TPE	SOFT-12-2BV
Large Adult Long (12L)	32-43 cm	SOFT-12L	SOFT-12L-1SC	SOFT-12L-1TPE	SOFT-12L-1MQ	SOFT-12L-1HP	SOFT-12L-2SC	SOFT-12L-2MQ	SOFT-12L-2TPE	SOFT-12L-2BV
Thigh (13)	40-55 cm	SOFT-13	SOFT-13-1SC	SOFT-13-1TPE	SOFT-13-1MQ	SOFT-13-1HP	SOFT-13-2SC	SOFT-13-2MQ	SOFT-13-2TPE	SOFT-13-2BV

\* For one male and one female locking fitting, change "MQ" to "MF" on the part numbers.

## Truly standardise your facility and stop worrying about tubes and connectors

- 1) Equip all blood pressure devices with FlexiPort fittings
- 2) Assign one cuff to patients upon admittance to the hospital
- 3) Have the unique ability to attach patients' cuffs directly to all blood pressure devices – room to room, floor to floor, and department to department

### Neonatal Soft Disposable Cuffs with Male Luer Slips

Size	Range	One-Tube (10 pack)	Two-Tube (10 pack)
NEO 1	3.3-5.6 cm	5082-101-1	5082-101-2
NEO 2	4.2-7.1 cm	5082-102-1	5082-102-2
NEO 3	5.4-9.1 cm	5082-103-1	5082-103-2
NEO 4	6.9-11.7 cm	5082-104-1	5082-104-2
NEO 5	8.9-15.0 cm	5082-105-1	5082-105-2
Multi-Pack	1 of each	5082-241-9	5082-241-10

For more information, contact your local Welch Allyn representative or visit us online at [www.welchallyn.com](http://www.welchallyn.com).

#### UK Customer Services

Welch Allyn UK Ltd, Cublington Road, Aston Abbots, Buckinghamshire HP22 4ND, England  
Tel: +44 (0) 207 365 6780 Fax: +44 (0) 207 365 9694

Welch Allyn South Africa (Pty) Ltd, 394 Surrey Avenue, Ferndale, Randburg, 2194, South Africa  
P O Box 809, Fountainebleau, South Africa 2032  
Tel: +27 (0) 11 777-7555 Fax: +27 (0) 11 777-7556

#### Welch Allyn Customer Services

Navan Business Park, Dublin Road, Navan, Co. Meath, Republic of Ireland  
Tel: +353 46 90 67790 Fax: +353 46 90 67754



Advancing Frontline Care™