



3.5 V LITHIUM ION RECHARGEABLE HANDLE



Twice the On-Time, Half the Weight. The Welch Allyn Power Handle...Smarter Technology.

Welch Allyn's smarter Lithium Ion technology provides twice the on-time with just half the weight of conventional 3.5 V handles. The Lithium Ion Power Handle's technology powers any of our 3.5 V instrument heads for 110 minutes—twice as long as the conventional handle. With an ergonomic design and weighing only 4.4 ounces, the handle is comfortable for any size hand. An indicator light warns you when the battery is running low, and your choice of charging modules makes recharging simple and convenient. Made of durable CYCOLOY®, the Lithium Ion Power Handle will provide years of reliability.

WelchAllyn®

Advancing Frontline Care™

Lithium Ion Technology Provides Twice the On-Time with Just Half the Weight

- Half the weight of conventional 3.5 V handles
- Smallest diameter of any full-size 3.5 V power handle
- Soft-grip rheostat intensity control section
- Textured handle
- Average 110 minutes of on-time between charging
- Low-battery indicator light
- Lithium Ion Battery Technology
 - No “memory” effect
 - Low battery indicator
 - No need to discharge between charging cycles
- Battery section built to last using CYCOLOY® resin

*CYCOLOY is a registered trademark of General Electric Company
- ETL listed to UL 2601 and CAN/CSA C22.2 No. 601.1 medical electrical equipment standards

(Handle shown at actual size)

Choice of Charging Modules:



Removable, Pull-out AC Charging Module

- Easily removable from handle
- Retractable AC prongs
- Can be attached or detached to handle during exams

71950 AC Charging Module only



Lithium Ion Handle Well Adapter

- Fits Welch Allyn 71110 desk charger
 - “Right sizing” top adjusts to different well heights
- 71980 Lithium Ion Well Adapter only (for use with 71110 Desk Charger and most Stand Console Chargers)

Ordering Information

71900 Lithium Ion Power Handle with AC Charging Module

71930 Lithium Ion Power Handle with Well Adapter

71910 Handle without AC Charging Module

71950 AC Charging Module only

71980 Lithium Ion Well Adapter only (for use with 71110 Desk Charger and most Stand Console Chargers)

71990 Lithium Ion Adapter only (for use with select Stand Console Chargers)

71960 Replacement battery

U.S. Patent No. 6,130,520

For more information about the Welch Allyn Lithium Ion Power Handle, contact Welch Allyn at 888.359.2394.

4341 State Street Road, PO Box 220, Skaneateles Falls, NY 13153-0220 USA
(p) 888.359.2394 or 315.685.2679 (f) 315.685.2174 www.welchallyn.com

© 2006 Welch Allyn SM2253 Rev C

WelchAllyn®

Advancing Frontline Care™